



National Institute for Occupational Safety and Health
National Personal Protective Technology Laboratory
Technology Evaluation Branch
Certification, Evaluation and Testing Section
P.O. Box 18070
Pittsburgh, PA 15236

TEST REPORT

Task Number: TN-19825

Manufacturer: Guangzhou Harley Commodity Company Limited

Prepared by: Nichole Petitta

Tests Conducted by: Nichole Petitta

Date: November 21, 2014

Respirator Tested: L-288 N95

Background Information

Respirator tested as per test request.

Tests Assigned

<u>Test Description</u>	<u>STP Number</u>
A. Exhalation Resistance Test	RCT-APR-STP-0003
B. Inhalation Resistance Test	TEB-APR-STP-0007
C. Sodium Chloride (NACL) N95 Test	TEB-APR-STP-0059

Overall Results

The respirator system tested met the requirements of all the procedures listed above.

Task Number: TN-19825

Reference No.: CFR 84.180

Test: Inhalation Resistance Test

STP No.: 7

Manufacturer: Guangzhou Harley Commodity Company Limited

Item Tested: L-288 N95

Comments:

Samples were tested on Brooks meter 000036.

Was all equipment verified to be in calibration throughout all testing?

Yes

No

Signature:

Nichole L. Petitta

Date: 11/4/2014

Engineering Technician

Task Number: TN-19825

Reference No.: CFR 84.181

Test: Sodium Chloride (NaCl) - N95

STP No.: 59

Manufacturer: Guangzhou Harley Commodity Company Limited

Item Tested: L-288 N95

Comments:

Samples were tested on TSI 8130 machine 000192 with samples 1-3 using timer 000214.

Was all equipment verified to be in calibration throughout all testing?

Yes No

Signature:

Nichole L. Petitta

Date: 11/21/2014

Engineering Technician